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
文書の許諾 Grason-Stadler, Inc.

Information in this Application Note is based on personal communications with, and the published research of, the following: Professor Field Rickards and colleagues in the Department of Otolaryngology, The University of Melbourne, Australia; Dr. Barbara Cone-Wesson, Department of Speech & Hearing Sciences, University of Arizona, Tucson; and Mr. John Parker, Interrupt Systems, Victoria, Australia.

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- Rickards, F.W.
Auditory steady-state evoked potentials in humans to amplitude modulated tones.
Ph.D. Thesis, August 1983.
Department of Otolaryngology, The University of Melbourne.
- Rickards, F.W. and Clark, G.M.
Steady-state evoked potentials to amplitude modulated tones.
Evoked Potentials II, edited by R.H. Nodar and C. Barber.
Butterworth, Boston, MA, 1984; 163-168
- Cohen, L.T.
Auditory steady-state evoked potentials:
Automated threshold estimation in awake and sleeping subjects.
Ph.D. Thesis, July 1990.
Department of Otolaryngology, The University of Melbourne.
- Cohen, L.T., Rickards, F.W., Clark, G.M.
A comparison of steady-state potentials to modulated tones in awake and sleeping humans.
Journal of the Acoustical Society of America, 1991; 90:2467-2479.
- Rance, G., Beer, D.E., Cone-Wesson, B., Shepard, R.K., King, A., Rickards, F.W., and Clark, G.M.
(1999) Clinical findings for a group of infants and young children with auditory neuropathy.
Ear and Hearing, 20 (3), 238-252.
- Tan, Lesley
Establishing and Comparing Thresholds in Babies using SSERA and ABR, Paper presented at the 13th National Conference of the Audiological Society of Australia, Sydney, 29 April to 1 May 1998.
- Burton M.J., Cohen L.T., Rickards F.W., McNally K.I., Clark, G.M.
Steady-state evoked potentials to amplitude modulated tones in the monkey.
Acta Otolaryngol (Stockh) 1992; 112:745-731.
- Rance, G., Rickards, F.W., Cohen, L.T., Burton M.J., Clark G.M.
Steady-State Evoked Potentials: A New Tool for the Accurate Assessment of Hearing in Cochlear Implant Candidates.
Advances in Oto-Rhino-Laryngology, Basel, Karger, 1993; 48:44-48.
- Rickards, F.W., Tan L.E., Cohen L.T., Wilson O.J., Drew J.H., and Clark G.M.
Auditory steady-state evoked potentials in newborns.
British Society of Audiology, 1994; 28: 327-337.
- Rance, G., Rickards, F.W., Cohen, L.T., De Vidi S., and Clark G.M.
Automated prediction of hearing thresholds in sleeping subjects using auditory steady-state evoked potentials.
Ear and Hearing, 1995; 16:499-507
- Rance, G., Dowell, R.C., Rickards, F.W., Beer, D.E., and Clark G.M.
Steady-state evoked potential thresholds in a group of children with absent click-evoked auditory brainstem response.
Ear and Hearing, 1998, 19:48-61.

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© 2001 Grason-Stadler Inc.