




# IMPEDANCE AUDIOMETER RS -H 1



# RS-H1

- Multifunctional type which is suitable for clinical and research corresponding to numerous tests
- Automatic measurements for a variety of test from tympanometry to reflex tests.
- The test results are stored in the internal memory and printed out in batches for each tested item, inspected ears.



1000 Hz probe tone, can do tympanometry that is suitable for newborns.

For both 226 Hz and 1000 Hz, the Tympano-gram can be compared with right and left or 5 data.

Equipped with LAN interface, also supports card reader and bar code reader

# Multifunctional type for clinical and research with enhanced function and operation.

## Characteristic

- The test can select either automatic start or manual start, and can execute the automatic measurement of both tympanometry and reflex test consecutively.
- User can change various condition settings such as sweep stop pressure of tympanometry, the start tone pressure level, level step of Reflex test.
- The test result is stored in the internal memory and is printed out by batch for each test item, test ears.
- In addition to the built-in servo pump, pressure setting of the reflex test is possible by using an external mechanical pump (accessory).

### Smaller size to save space

The width is 15.7% (about 6.5 cm) and the depth is 24.9% (about 9.8 cm) smaller than our previous product (RS - 22)\*.

Our same previous product \*



### Easy operation with touch panel

Some actions such settings and ID input can be operated with the such as touch panel, making it easier to use.



### Easier to use and improved ear probes

The transition from the handheld to the headset is smooth.



### Mechanical pump

#### FB-22



Mechanical pump that can measure pressure with reflex test is included standardly

### Easy to use operation panel

Equipped with a probe holder in the main body



### Equipped with LAN interface


Test data can be transferred to external equipment such as a computer via LAN interface. (Please contact our sales department for use)

## Tympanometry


- It is equipped with 1000 Hz probe tone. It is possible to do Tympanometry which is suitable for newborn babies.
- The left and right test results can be displayed in color on the same screen, and 5 patterns can be displayed on the same screen.
- Automatically measure and display ear canal volume, static compliance, and peak pressure. (Probe tone 226 Hz).
- Automatically measure and display ear canal volume, static compliance, and peak pressure. Static acoustic admittance can be displayed by selecting from "Positive Tail" or "Negative Tail". (Probe tone 1000 Hz)
- "Gradient" showing the steepness of the peak can also be displayed.
- Tympanometry "Relative display (MCT)", "Absolute display (MPT)" can be switched.

### Tympanogram

Probe tone 226 Hz




Probe tone 1000 Hz




### Printout example

Left and right test result simultaneously



5 pattern simultaneously






## Reflex test


- Five modes (AUTO, MANUAL, ADD, DECAY, LATENCY) are carried.
- Easy threshold measurement in automatic mode (AUTO), measurement with arbitrary level, frequency, timing by reflex test (MANUAL) is possible.
- Easy reflex test (DECAY) is possible with 10 seconds stimulation.
- In addition to "real wave punishment display", "simple display" can also be selected. (Reflex test <Auto> only) Measurement of latency time (Reflex check (Latency)) has become more convenient.
- It is possible to independently set the cursors of the same side waveform and the opposite side waveform. And it is possible to display the latency time of each side

### Test mode: AUTO



For the set frequency, measure the same side / opposite side automatically by changing the stimulus tone.


### Test mode: MANUAL




Measure with the set frequency, level, timing

## Printout examples


### Reflex test (Auto)



### Reflex test (Manual)




## Test mode: ADD




Measure by adding at the set frequency, level, timing. Select from adding time 5, 10, 20.

## Test mode: DECAY



Measure with posting stimulus tone for 10 seconds at the set frequency and level. Measurement time 12.5 seconds

## Test mode: LATENCY




Measurement is performed by adding and averaging at the set frequency, level and timing. Select adding time from 1, 5, 10, 20 times. Time measurement at an arbitrary place is possible

## Cooperation with the electronic medical record system

User can transfer test results directly by using optional data linkage system.

## Patient Information


Patient information can be recalled from the electronic medical record system by ID number. The ID number can be entered with an optional card reader or bar code reader.




Selection of the audit's name

## Image output is also possible.

Images can be transferred directly to the computer without using the data linkage system.



## Electronic medical record system



Please contact our sales department for use.

# IMPEDANCE AUDIOMETER RS-H1

## Specification

Type according to IEC60645—5: Type 1—Diagnostic／Clinical

### Test function

#### Tympanometry

Probe tone: 226 Hz, 85 dB SPL, 1000 Hz 75 dB SPL (According to IEC 60318 - 5: 2 cm<sup>3</sup> coupler for the measurement of hearing aids and earphone coupled to the ear by means of ear inserts).

Display mode: Left and right overwriting or 5 data overwriting.

Equivalent volume range: When probe tone 226 Hz is selected.

0.10 ~ 9.99 mL (display range).

0.20 ~ 9.00 mL (error guarantee range).

When probe tone 1000 Hz is selected.

0.22 (0.05 mL) ~ 35.00 (8.00 mL) mmho (display range).

0.44 (0.10 mL) ~ 26.50 (6.00 mL) mmho (guaranteed error range).

Range: mL / div 0.5, 1.0, 1.5, 2.0, 3.0 Automatic selection functions are available.

mmho / div 0.5; 1.0; 2.0; 4.0; 10.0 Automatic selection functions are available

Pressure range: + 200 ~ -600 daPa; + 200 ~ -400 daPa; + 200 ~ -200 daPa; + 200 ~ automatic stop.

(Automatic stop of conditions: - 200 daPa equivalent volume value in the following stop when down to about 1/3 of the static compliance. Maximum -400 daPa).

The screen format has two types of pressure range + 200 ~ -600 daPa and + 200 ~ -400 daPa

Pressure group acceleration: 200 daPa / sec, 50 daPa / sec.

Display: When 226 Hz is selected.

Tympanogram can be overwritten left right or overwritten with 5 data.

Numerical values are peak pressure, S. C, PVT, Gradient, Probe tone frequency.

Relative display (Meatus compensated Tympanometry).

Absolute display (Meatus compensated Tympanometry).

When 1000 Hz is selected.

Numerical values are peak pressure, SAA (p / n), ECV, Gradient, probe tone frequency.

Relative display (Meatus compensated Tympanometry).

Absolute display (Meatus compensated Tympanometry).

#### Reflex test

Probe tone: 226 Hz, 85 dB SPL

Test mode: The following five types

AUTO: Measured with the level of stimulus tone automatically changed and continuously on the same side / opposite side for the set frequency.

Threshold can be set

Can set the stop level of stimulus tone

MANUAL Measure with arbitrary stimulus tone frequency, level, timing

Automatic (Measured with the level automatically changed for the set frequency)

Threshold can be set

ADD Measure by adding at any stimulus tone frequency and level

Addition time can be selected from 5; 10; 20 times

DECAY Measure by posting stimulus tone for 10 seconds at arbitrary 1 stimulation tone frequency and level

Measurement time 12.5 seconds

LATENCY Measure by adding averagely at arbitrary 1 stimulation tone frequency and level

The addition time can be selected from 1; 5; 10; 20 times.

Has cursor function

Waveform display: Can switch from actual waveform to simplified display (AUTO only).

Sensitivity: 0.025 (AUTO only), 0.05, 0.1, 0.2 mL / div

Pressure range: + 200 ~ -600 daPa

Measurement pressure: atmospheric pressure, tympanogram peak pressure, manually set pressure.

Setting change of mechanical pump / servo pump is available.

Stimulus tone: As shown in the table above.

\* Specifications are subject to change without notice for improvement.

\* Data such as test screens and printout examples posted on this page are prepared for catalogs.

Frequency (Hz)	Pure tone						Noise			
	250	500	1000	2000	4000	8000	Wide Noise	Low Noise	High Noise	
Synchronism	(dB SPL)	—	50~110	50~110	50~110	50~110	—	50~90	50~90	50~90
	(dB HL)	—	50~105	50~110	50~105	50~105	—	—	—	—
Opposite-side	(dB SPL)	—	—	—	—	—	—	50~120	50~120	50~120
	(dB HL)	50~100	50~120	50~120	50~120	50~120	50~90	—	—	—

However, at the time of LATENCY test, the same side 500 Hz: 50 ~ 100 dB SPL, 50 ~ 90 dB HL

shall be set. Or, the noise output of the same side is not performed

Wide noise: 500 ~ 4000 Hz noise

Low noise: 500 ~ 2500 Hz noise

High noise: 2500 ~ 4000 Hz noise

Synchronization, opposite coupler and HK reference level are according to IEC 60645-5.

Minimum level step: 5 dB

#### Reflex (AUTO)

##### Stimulus tone

Frequency On both the frequency side and the opposite side, there are four frequencies of 500,

1000, 2000, 4000 Hz

ON / OFF setting is possible for each frequency

Level step 5 dB (7 level), 10 dB (4 level) can be selected

Start level 70 dB, 80 dB, optional

The threshold determination function 0.025 mL; 0.05 mL; OFF can be selected

Expanded display : Expand and display the waveform of the selected frequency.

Stop stimulation tone: Can set the stop level of stimulation tone.

#### Reflex (Manual)

##### Stimulus tone

Frequency: Same side or opposite side can be selected arbitrarily

Tone posting 1.5 seconds, manual (posted while pushing stimulus notice bulletin), can be selected

automatically (rising automatically in steps of 5 dB from the set start level)

Continuous measurement result (overwrite), single (for one screen only)

Stimulus tone: single only when automatic

Sweep speed 2.5; 5.0; 10.0; 20.0 seconds / screen selectable

Threshold judgment function 0.025 mL, 0.05 mL, OFF selectable

However, stimulation tone: valid only when automatic

#### Reflex (ADD)

Stimulus tone posting 0.5; 1.0; 1.5; 2.0 seconds selectable

Adding times 1; 5; 10; 20 times selectable

Sweep speed 2.5; 5.0; 10.0; 20.0 seconds / screen selectable

#### Reflex (LATENCY)

Stimuli tone posting 0.5; 1.0; 1.5; 2.0 seconds selectable

Adding times 1; 5; 10; 20 times selectable

Sweep speed 2.5; 5.0; 10.0 seconds / screen selectable

Cursor function

#### Tympanometry & flex test

Continuous test of Tympanometry 226 Hz (left and right over writing) and

Reflex Auto test

#### Digital section

Interface: RS-232-C: Data output, card reader / bar code reader

LAN: Data input / output

External switch: Start / stop measurement by switch

Waveform output: Electrical output

Indicator: 7 inch wide color LCD, 800 x 480 dots, resistive touch panel Built-in printer: 4 inch

thermal printer; 832 dots / line, 8 dots / mm

#### Power supply part

Classification of protection against electric shock by type: Class I device

Classification of mounting parts depending on level of protection against electric shock: B type mounting part

Power: AC100V 50/60 Hz 60VA

#### Environmental condition

Use environmental conditions: Temperature 15 ~ 35°C,

Relative humidity 30 ~ 90% (without any condensation)

Pressure 700 ~ 1020 hPa (altitude 3000m or less)

Transport and storage environmental conditions: Temperature -10 ~ 50°C,

Relative humidity 10 ~ 90% (without any condensation)

Pressure 500 ~ 1050 hPa

#### Size and weight

Body: About 350 (width) x 295 (depth) x 290 (height) mm about 6.0 kg

#### Accessories

Ear Probe	ET-05	1	Earplug set	RS-M1-S10	1
Earplug adapter for ET - 05	ET-05-S11	2	• Earplug for the RS		Each 3
Air conduction handset	AD-06BF2	1	• Big, Medium, Small		Each 3
Air conduction adapter assembly	EB-01A-017	1	• Earplug Φ 7~14		Each 3
Headband adapter for ET-05	ET-05-S12	1	• Earplug Φ 6~9		
Head band assembly for ET 05	ET-05-S13	1	• Cleaning Wire		
Fuse	0218002.MXP	2	Recording paper axis	RS-32-026	1
Mechanical pump	FB-22	1	Thermal recording paper	AA-M1-025	1
			3P power cord with ground	TP-42	2
				AA-38-222	1

Specifications subject to change without notice



RION CO., LTD.

Head Office / Sales Office

〒185-8533

3-20-41, Higashimotomachi, Kokubunji-shi, Tokyo

TEL: 042-359-7880 FAX: 042-359-7441

<http://www.rion.co.jp/english>